



ELECTROPHORESIS AUTOMATED SYSTEM

ADALYA IFE



ELECTROPHORESIS AUTOMATED SYSTEM



VERSATILE

Allows to perform the most common techniques: Serum proteins, Lipoproteins, Urinary Proteins, Immunofixation.

COMPLETE AUTOMATION

AdalyaIFE automatically arranges the loading of all the reagents, to standardize the delicate electrophoretic process.




Serum Proteins



Haemoglobins



Lipoproteins



Integrated webcam for reporting of results



Automated Ag/Ab reaction contact

PRODUCTIVITY

Serum-Proteins : 40 samples

Lipoproteins: 40 samples

Haemoglobins: 40 samples

Immunofixation: 5 samples

SOFTWARE FEATURES

User interface: monitor, keyboard and mouse

Minima recognition: automatic

Managing techniques: adjustable techniques in times, migration voltage/current

Operating system:



ADALYA IFE



SELEO
Engineering

Your partner in Protein Electrophoresis