

# MiniCute Imager

Small, Light, and Easy

- **Switchable filters:**  
Switchable UV and blue light filters
- **Compact:**  
Smaller dimension, save your working area
- **Convenient:**  
Wi-Fi compatible, easy to download and print the image
- **Easy-to-use:**  
No complex interface, just operate the digital camera
- **Open system:**  
High compatibility to other brand transilluminators



MiniCute Imager is an easy gel image system. Compatible with MaestroGen UltraSlim® LED, UltraBright LED and UltraBright UV Transilluminators. Embedded with Canon digital camera, 100% dark room for photography.



## Specification:

Product	MiniCute Imager
Cat No.	MC-02
Outer dimensions	WxDxH 275x238x250mm
Filter	Amber filter (for blue light) and UV-cut filter
Weight	2 kg
Digital camera	Canon EOS M-series Type (1) Effective pixel: 24.7 Mega pixels (2) LCD Size: 3" TFT Color (3) Lens: 5x lens, F1.8 (Wide) – F8.0(Telephoto) (4) Mode: Program (P), Aperture Priority (A), Shutter Priority (S), Manual (M) (5) ISO: Auto/80/100/125/160/200~25600 (6) Shutter Speed: 30~1/4000 second (7) Video Output: NTSC / PAL selectable (8) Interface: USB (v.2.0), Wi-Fi compatible (9) Storage Media: SD/SDHC memory card (10) Storage Format: JPEG (Exif Ver 2.30)/RAW (11) Weight: 320g (12) Dimensions: 108.2x67.1x35.1mm

## Related product:

Product	UltraSlim® LED Illuminator	UltraBright LED Transilluminator	UltraBright UV Transilluminator
Cat No.	SLB-01 /SLB-01W	LB-16	MLB-21/MB-21 MLB-16/MB-16
Description	viewing stained gels on the laboratory bench, or within gel documentation system		