



# Hematology Analyzer

Model: Nano 3



## 3 modes

Venous whole blood.  
Capillary whole blood.  
Pre - diluted blood.

- 3-part differentiation.
- 9 uL sampling volume.
- Up to 60 samples per hour.
- 21 parameters.



Arasaat Al-Hindia Dist. 929 St.21, Villa 32  
Baghdad-IRAQ

0783 330 0730 البصرة 0771 608 6395 بغداد  
0751 815 6446 الموصل 0751 815 6446 أربيل  
0770 562 46 69 العمارة 0770 681 5747 سلیمانیة  
0790 238 8611 الأنبار 0770 135 0023 كركوك  
0774 024 2223 واسط

شركة أفلو للأجهزة الطبية و المختبرية

afloltd

@afloltd

www.afloltd.com



شركة أفلو  
للأجهزة الطبية و المختبرية



## Operation Principles

Electrical impedance method (WBC/RBC/PLT) and cyanide free (HGB) method

Pulse baseline tracking technology, waveform discrimination technology, digital sheath flow technology and intelligent

floating threshold technology

Dual-channel is used for counting

## Measured Parameters

WBC, Lym%, Mid#, Gran#, Lym%, Mid%, Gran%, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-CV, RDW-SD, PLT, PDW, MPV, PCT, P-LCC, P-LCR, 3 histograms, alarm for abnormal erythrocytes, leukocyte, and platelet

## Throughput

Up to 60 samples per hour

## Hole blood sample

Minimum volume 9 ul

**Digital Bank** : Male, Female, New born Control blood.

**Modes** : Verous whole blood , anticoagulant periphal whole blood and pre-diluted blood

**Analysis time** : Results in 40 Seconds

**Display** : 11" Touch Screen TFT display

## Performance

Carryover Repeatability Linearity

WBC  $\leq 0.5\% \leq 2.0\%$  (7-15 $\times 10^9$ /L) 0.00-300 $\times 10^9$ /L

RBC  $\leq 0.5\% \leq 1.5\%$  (3.5-6.0 $\times 10^{12}$ /L) 0.00-8.5 $\times 10^{12}$ /L

HGB  $\leq 0.5\% \leq 1.5\%$  (110-180g/L) 0-250g/L

PLT  $\leq 1.0\% \leq 4.0\%$  (150-500 $\times 10^9$ /L) 0.00-3000 $\times 10^9$ /L

## Data Storage

50,000 sample results, including parameter, histogram and patient information

**Languages** : English ,French , Spanish, German, Italian, Chinese

## Printer

Internal Printer with USB Ports and LIS connection

**Operating Temp:** 9C to 35C

**External Power Supply** 100-240Vac , Power 300VA, Frequency :50/60Hz

**Dimensions:** 15.5"(L)  $\times$  14"(W)  $\times$  18.5"(H)

**Weight:** 23.8kg

**Reagents** : Diluent 3-PAC

Lyse: ISOFLEX

Cleaner: ACT-CLEAN

